

Indirect ELISA

All steps are carried out at room temperature unless otherwise indicated. Recipes for all solutions are listed at the end of the protocol.

Antigen Coating

1. Dilute purified antigens to a final concentration of 0.2 µg/ml in antigen-coating buffer and add 100µl of diluted antigen to each well of a 96-well ELISA plate.
2. Carefully cover the plate with adhesive plastic and incubate at 4°C overnight.

Blocking

3. Empty the wells of antigen-coating buffer and wash twice with 200 µl PBST buffer for 5 minutes each time.
4. Add 200 µl blocking buffer per well to block residual protein-binding sites. Cover the plate with adhesive plastic and incubate for 1-2 hour(s) at 37°C.

Antibody incubation

5. Dilute your primary antibody of choice with blocking buffer in a series e.g. 1:500, 1:1000, 1:2000, 1:4000 and so on, empty the wells of blocking buffer and then add 100 µl of each dilution per well. Repeat in duplicate, or triplicate, for accuracy. Cover the plate with adhesive plastic and incubate for 1 hour at 37°C.
6. Empty the wells and wash 3 times with 200 µl PBST buffer for 5 minutes each time.
7. Dilute the HRP-conjugated secondary antibody with blocking buffer at an optimal concentration (a dilution factor within 1:10,000-1:100,000 is recommended) and add 100 µl of secondary antibody solution to each well. Cover the plate with adhesive plastic and incubate for 1 hour at 37°C.
8. Empty the wells and wash 3 times with 200 µl PBST buffer for 5 min each time.

Antibody incubation

9. Add 100 µl TMB substrate (mix equal volumes of TMB buffer A and buffer B) to each well with a multichannel pipette. Color development should peak after 15 minutes, at which time it should be stopped by adding 100µl of 2M H₂SO₄ per well. Read absorbance at 450nm.

Solutions

Blocking buffer (100 ml)

5% non-fat dry milk 5g, add PBST buffer to 100ml

Antigen-coating buffer (1000ml)

100mM NaHCO₃ 8.4g, adjust pH to 9.6, add ddH₂O to 1000ml

PBST buffer (1000ml)

10mM Na ₂ HPO ₄	1.42g
1.8mM NaH ₂ PO ₄	0.22g
140mM NaCl	8.19g
0.2% Tween 20	2ml
Adjust pH to 7.4	
Add ddH ₂ O to 1000 ml	

TMB buffer A (500ml)

NaAc•3H ₂ O	13.6g
Citric Acid	1.6g
30% H ₂ O ₂	0.3ml
Add ddH ₂ O to 500 ml	

TMB buffer B (500ml)

TMB (dissolved in 3ml DMSO)	0.15g
EDTA-2Na	0.2g
Citric Acid	0.95g
Glycerol	50ml
Add ddH ₂ O to 500 ml	

For research purposes only ! Not for therapeutic or diagnostic purposes in humans or animals !